

Hue Saturation Brightness workflow tutorial

Step by step instructions

Instructions are correct as of 21/05/18.

Links to websites are subject to change and not guaranteed. It is assumed that users will seek the appropriate rights, licenses or permissions to download, distribute and use them according to the relevant copyright laws from the respective owners. We do not own or endorse any of the software mentioned, and this tutorial is strictly for the purpose of illustrating the HSC workflow only.

This tutorial assumes that you have a competent 64-bit computer running **Windows 10**, and have the necessary administrator rights. For the purposes of this tutorial, it also assumes that the computer has a fully licensed version of **Imaris (version 8.4.2)**, **FlowJo VX (version 10.4.1)** and **Microsoft Excel (2010)**. In many cases, older or newer versions of the software should have the same functionality as described, but may not necessarily be accessed in the same manner described.

For clarity, in this tutorial, full filepaths are always given assuming a user account named “Lenovo” (**C:\Users\Lenovo**) on the C drive. Users should substitute the file paths accordingly.


Table of Contents

1	Getting Started	3
1.1	Initial Set Up of software environment for Windows 10, FIJI, WinPython 3.5.1.3 and WinPython 2.7.10.2.....	3
1.2	Initial Set Up of software environment for HSC scripts	4
1.3	Initial Set up of software environment for running Python on Imaris	5
1.4	Loading HSC LUTs into FIJI and Imaris.....	7
1.5	Test Data	7
2	Using the Hue-Saturation-Brightness surface creation (HSC) workflow.....	10
2.1	Rename TIF files to standard format	10
2.2	Reduce number of channels to 3 by merging channels (Optional).....	11
2.3	Generate Hue, Saturation and Brightness channels	11
2.4	Surface creation using the Brightness channel in Imaris and statistics extraction. ..	14
2.5	Saving out Surface Statistics.....	19
2.6	Aggregate statistics into single CSV file and Gating in FlowJo.....	20
2.7	Backgating in Imaris (Version 8.4.2).....	27
3	Spectral Compensation (optional).....	30
3.1	Compensation using FlowJo	30
3.2	Correction of original images.....	37

1 Getting Started

1.1 Initial Set Up of software environment for Windows 10, FIJI, WinPython 3.5.1.3 and WinPython 2.7.10.2.

1. Press the **Windows button** and immediately type “Show hidden files and folders”. Click on the corresponding search result to display the **File Explorer Options** window. Ensure that the option for **Hide extensions for known file types** is left unchecked. Click **OK** to save the settings.



Unhiding extensions for known file types

▲Critical step

This ensures that you are always able to recognize the file type (e.g. .py, .txt, .tif etc.) of the file that you are working on.

2. On the desktop (*C:\Users\Lenovo\Desktop*), create a folder called *Software*.
3. Create a subfolder called *FIJI* in the folder *Software*.
4. Download **FIJI is Just ImageJ (FIJI)** from the source website (<https://fiji.sc/>). Unzip the contents of the FIJI download into this folder.
5. Find the file *C:\Users\Lenovo\Desktop\Software\FIJI\fiji-win64\Fiji.app\ImageJ-win64.exe*. Select the file, right click, and choose **Send to > Desktop (create shortcut)**. This creates a shortcut on the desktop that allows you to launch **FIJI** conveniently.
6. Create another subfolder called *WinPython35* in the folder *Software*.
7. Download the 64 bit **WinPython (version 3.5.1.3)** from the source website (https://sourceforge.net/projects/winpython/files/WinPython_3.5/3.5.1.3/). Do not

download the “Zero” versions. Run the installer and follow the instructions. When prompted for the install location, browse to the *C:\Users\Lenovo\Desktop\Software\WinPython35* folder.


► Troubleshooting

If you are running **Windows 7** or **8** instead of **Windows 10**, note that you might also need to install **Microsoft Visual C++** packages in order to get **WinPython** to work. See <https://github.com/winpython/winpython/wiki/Installation> for details (e.g. **Microsoft Visual C++ Redistributable for Visual Studio 2015** is required for **WinPython 3.5.1.3**. Use your search engine to find the required downloads from the official websites).

8. Create another subfolder called *WinPython27* in the folder *Software*.
9. Download the 64 bit **WinPython (version 2.7.10.2)** from the source website (https://sourceforge.net/projects/winpython/files/WinPython_2.7/2.7.10.2/). Run the installer and follow the instructions. When prompted for the install location, browse to the *C:\Users\Lenovo\Desktop\Software\WinPython27* folder.

1.2 Initial Set Up of software environment for HSC scripts

10. Create a subfolder called *HSC* in the folder *Software*.
11. Download the HSC scripts and look up tables from the supplementary information section in the manuscript. Unzip the contents into this folder. You should now have two new folders named *Data* and *Scripts* directly in the *HSC* folder.
12. In **Windows**, locate the first script *C:\Users\Lenovo\Desktop\Software\HSC\Scripts\P01_Rename.tif* files into standard **format.py**. To let **Windows** associate .py files with **WinPython 3.5.1.3**, right-click the file and select **Open with > Look for another app on this PC**, and browse to *C:\Users\Lenovo\Desktop\Software\WinPython35\WinPython-64bit-3.5.1.3\python-3.5.1.amd64\python.exe* to run the application. Ensure that the option **Always use this app to open .py files** is checked. Click **OK**.



Initial steps for associating .py files with python.exe


Since we are not yet doing anything, press **OK** when the program message appears, but then click **Cancel** when the program asks for the folder to process. Some error messages should pop up in a new command window, but you can safely ignore and close them.

□ Explanatory note

You should now be able to directly run .py scripts by double clicking on the file.

1.3 Initial Set up of software environment for running Python on Imaris

13. Open the **Imaris** software. To set up **Python 2.7** for Imaris XT, go to **File > Preferences > CustomTools**, click **Browse** under **Python 2.7 Application** and browse to *C:/Users/Lenovo/Desktop/Software/WinPython27/WinPython-64bit-2.7.10.2/python-2.7.10.amd64/python.exe* as illustrated in the screenshot below.




Setting up WinPython 2.7.10.2 to work with **Imaris XT**

14. Under the **XTension Folders** tab, click **Add...**, and browse to *C:\Users\Lenovo\Desktop\Software\HSC\Scripts\XT*. Click **Select Folder**. Click **OK** to close the settings window. Close the **Imaris** window and restart the software. Now, you should be able to see **Image Processing > XTIsolateSurfacesWithCsvList** appearing as an option.

► Troubleshooting

Click on **Image Processing > XTIsolateSurfacesWithCsvList** to test if the script is working. You should see a command window with a message asking you to create some surfaces in the Surpass scene first. However, if you do not see this screen, but instead get a command window that flashes briefly and then immediately disappears, the XT is not working. Try changing the filepath in **Step 18** to *C:/Users/Lenovo/Desktop/Software/WinPython27/WinPython-64bit-2.7.10.2/python.exe* instead. Contact your software manufacturer if this still does not work.



Pop up command window showing WinPython script working with **Imaris XT**

1.4 Loading HSC LUTs into FIJI and Imaris

In order to visualize the hue of surfaces and cells, a customized look up table (LUT) for the hue channel has been created for **FIJI** and **Imaris**.

15. Find the file *C:\Users\Lenovo\Desktop\Software\HSC\Scripts\HSB_Hue.lut*, and copy and paste the file into *C:\Users\Lenovo\Desktop\Software\FIJI\fiji-win64\FIJI.app\luts*. **FIJI** will automatically detect the LUT the next time it starts up.
16. Find the file *C:\Users\Lenovo\Desktop\Software\HSC\Scripts\HSB.pal*. Go to the folder where **Imaris** is installed, (e.g. *C:\Program Files\Bitplane\Imaris x64 8.4.2\colorTables*) and copy and paste the file in. **Imaris** will automatically detect the LUT the next time it starts up.

1.5 Test Data

17. Create a new folder called *Data* on the desktop (*C:\Users\Lenovo\Desktop\Data*). Copy and paste the raw data file (**Olympus** format) *C:\Users\Lenovo\Desktop\Software\HSC\Data\Testdata.oib* into this new folder.
18. Start **FIJI**, and go to **File > Open**. Select the filepath for *C:\Users\Lenovo\Desktop\Data\Testdata.oib* and click **OK**. Use the default import settings and click **OK**.

□Explanatory note

By default, **FIJI** uses the **Bio-Formats Importer** (under **Plugins > Bio-Formats > Bio-Formats Importer**) to read in data from most major microscope manufacturers.

□Explanatory note


With an image file open and selected, by mousing over the image, the **FIJI** main window will display the specific intensity value of the pixel where the mouse is pointing to, as well as other useful information about the image. An image is simply a grid of numbers (representing intensities), where the computer uses a look-up table (LUT) to map each number to a specified color to display on screen.

19. Go to **Image > Type** and select **32bit**.

□Explanatory note

This converts the images from 16 bit to 32 bit for compatibility with some scripts.

20. Go to **File > Save As > Image Sequence...** and the following dialog box should pop up.



Save Image Sequence Dialog box in **FIJI**

21. Select the format as “**TIFF**” and key in “**Testdata_c**” (i.e. without the literals/ inverted commas) under the **Name** category and click **OK**.

□ Explanatory note

For multidimensional data (i.e. hyperstack), **FIJI** typically names them in the following format: e.g. **XXX_t0001_z0005_c0009.tif**, and no further adjustments are normally necessary. However, in the case of a single dimensional data, its default behaviour is to only include the numbering after the prefix, e.g. **XXX0006.tif**. To disambiguate whether the number represents channel number (c), z slice position (z), or time index (t), we had to manually include “**_c**” in the filenames to inform the subsequent script that they represent channel numbers in this case.

▲ Critical step

It is important that the “**c**” be keyed in lower case and not as the capital letter “**C**”.

22. Create a subfolder **TIFF** under the **Data** folder, and save the TIFF images there. 4 TIFF images should be now saved under the **TIFF** folder, labelled **Testdata_c0000.tif** to **Testdata_c0003.tif**.

23. Only the files **Testdata_c0000.tif** to **Testdata_c0002.tif** contain fluorescent microscopy data. Delete **Testdata_c0003.tif**.

□ **Explanatory note**

Image Channel Information:

Testdata_c0000: DAPI

Testdata_c0001: CFSE


Testdata_c0002: Propidium Iodide

Testdata_c0003: Transmitted light (to be deleted)

2 Using the Hue-Saturation-Brightness surface creation (HSC) workflow

2.1 Rename TIF files to standard format

24. Go to *C:\Users\Lenovo\Desktop\Software\HSC\Scripts*, run the script named *P01_Rename tif files into standard format.py*. The script should start, and the following dialog box will appear:




The image shows a graphical user interface for the P01 script. It has a 'Directory:' field with a 'Select directory' button. Below it is a 'Match string:' field containing the regular expression `(t)(\d*)(_c)(\d*)(_xy)(_z)(\d*)(\.registered)(\.tif)`. There are several input fields for parameters: 'Prefix:' with value 'q', 'ttt:' with value '2', 'ccc:' with value '4', 'xxx:' with value '-1', 'yyy:' with value '-1', and 'zzz:' with value '7'. At the bottom are 'Done' and 'Cancel' buttons.

Graphical User Interface for P01

For more information regarding the parameters of each individual script, please refer to Supplementary Note 3 where a detailed description of each parameter is written.

25. Click **Select directory**, and select the file path for *C:\Users\Lenovo\Desktop\Data\TIFF*. The script should complete running, and the following dialog box will appear:



```
C:\Program Files\WinPython-64bit-3.5.3.1Qt5\python-3.5.3.amd64\python.exe
Directory chosen: C:\Users\Lenovo\Desktop\Data\TIFF
...Started processing. Please wait until program says it is done.
New folder called 'P01_Renamed TIFFs' created
...Started duplicating files. Please wait.
P01_Renamed TIFFs is not a file.
Copying "Testdata0000.tif"
Copying "Testdata0001.tif"
Copying "Testdata0002.tif"
Copying "Testdata0003.tif"
Working in new folder: C:\Users\Lenovo\Desktop\Data\TIFF\P01_Renamed TIFFs
P01_Rename TIFF files_log.txt could not be matched!
Renamed 'Testdata0000.tif'
Renamed 'Testdata0001.tif'
Renamed 'Testdata0002.tif'
Renamed 'Testdata0003.tif'
...Done. Press 'Enter' to exit.
```

Command window after P01 has completed running

26. Close the command window.
27. A new folder *P01_Renamed TIFFs* is now saved as a subfolder under the earlier selected path. A log file *P01_Rename TIFF files_log.txt* and three renamed TIFF files in the format *q_t0000_c00_x000_y000_z0000.tif* to *q_t0000_c02_x000_y000_z0000.tif* are found in this folder.

□ Explanatory note


The script renames the TIFF files it finds in the folder that matches a set of predefined nomenclature patterns. In this case, the script looked for the pattern “*XXX_c0000.tif*” and recognized that the numbers after “*_c*” represent the channel numbers. The TIFF files need to be renamed into a specific convention so that subsequent scripts can

correctly identify the position of each image. The script is designed to read in multiple types of filenames from **FIJI** and **Imaris**, and can also be modified to rename TIFF files with custom nomenclature from other sources. Refer to the manual for details.

2.2 Reduce number of channels to 3 by merging channels (Optional)

In our current test dataset, there is no need to merge channels, since we only have 3 channels - each be designated as one of the three RGB input channels. However, if you are using an image with 4 or more channels, you may need to merge channels. We have provided a script that calculates the maximum intensity projection of the input channels. We provide the instructions in this section as part of the full HSC workflow, but in this tutorial, we will not be using them.

28. Run the script **P02_Calculate max intensity channel.py**.



Graphical User Interface for P02

29. Key in the appropriate input channel numbers for merging (as a list) in the parameter **inChNumbers** and designate the output channel name with **outChNumber**.

▲Critical step

The number of each channel can be separated using a space or a comma.

30. Input the minimum values of each channel in the parameter **minValues** and the maximum values of each channel in the parameter **maxValues**. This allows each channel to be scaled according to its dynamic range. The user can decide on the values to input manually or they can obtain these values by opening the images in **FIJI**.

Pressing **Ctrl+Shift+C** in **FIJI** opens the brightness and contrast panel, and the program automatically scales each channel to its minimum and maximum values when the **Auto** button is pressed.

31. Click **Select directory** and choose the folder that contains the TIFF files (usually **P01**) and press **Done**.

▲Critical step

New files with the channel number **outChNumber** will be created in the same folder as the input path. Ensure that they do not have the same channel number as any other image files within the folder as they will be overwritten without warning.

2.3 Generate Hue, Saturation and Brightness channels

32. Run the script `P03_Generate HSB channels.py`. The following screenshot should show up:

Parameter	Value
inChNumberA :	02
inChNumberB :	01
inChNumberC :	00
sfA :	1.0
sfB :	1.0
sfC :	1.0
minA :	50
maxA :	3500
minB :	220
maxB :	3000
minC :	150
maxC :	3500
outChNumberH :	67
outChNumberS :	66
outChNumberB :	65

Graphical User Interface for P03

33. Key in the following values for the input channel numbers:

```
inChNumberA = "02" #'Red'  
inChNumberB = "01" #'Green'  
inChNumberC = "00" #'Blue'
```

□ **Explanatory note**

The parameters for A, B and C are always for red, green and blue respectively. Thus, in this case, channel “02” (Propidium Iodide) is designated red, channel “01” (CFSE) is green, and channel “00” (DAPI) is blue. The assignment of colors is purely cosmetic, but we used here an assignment that matched the actual dye colors.

34. Set the scale factors for the input channels using the following values:

```
sfA = 1.0  
sfB = 1.0  
sfC = 1.0
```


35. Set the minimum and maximum values for the input channels using the following values:

```
minA = 50  
maxA = 3500  
minB = 220  
maxB = 3000  
minC = 150  
maxC = 3500
```

36. Keep the output channel numbers as follows, and save the script before closing it.

```
outChNumberH = "67" #Hue  
outChNumberS = "66" #Saturation  
outChNumberB = "65" #Brightness
```

37. Click **Select directory** and choose the file path for *C:\Users\Lenovo\Desktop\Data\TIFF\P01_Renamed TIFFs*.
38. The following command window should appear. Close the command window.



```

C:\Program Files\WinPython-64bit-3.5.3.1Qt5\python-3.5.3...
Directory chosen: C:/Users/Tan Yingrou/Desktop/Data/TIFF/P01_Renamed TIFFs
...Started processing. Please wait until program says it is done.
New folder called 'P04_HSB' created
"P01_Rename TIF files_log.txt" was ignored.
"P02_Calculated Max Intensity Channel" was ignored.
"P04_HSB" was ignored.

Processing q_t0000_c02_x000_y000_z0000.tif, q_t0000_c01_x000_y000_z0000.tif, q_t
0000_c00_x000_y000_z0000.tif...
RGB to HSB ok
q_t0000_c67_x000_y000_z0000.tif was saved...
q_t0000_c66_x000_y000_z0000.tif was saved...
q_t0000_c65_x000_y000_z0000.tif was saved...

...Done.
Enter to exit.

```

Output command window from script P03

39. A new subfolder *P03_HSB* containing the log file and TIFF files of the **Hue** (“c67”), **Saturation** (“c66”) and **Brightness** (“c65”) channels is now saved under the file path *C:\Users\Lenovo\Desktop\Data\TIFF\P01_Renamed TIFFs\P03_HSB*.

□Explanatory note

For a two-channel image, for instance a GFP and YFP image, the ‘blue’ channel (**Channel C**) uses the YFP channel as a placeholder, and large numbers can be used as the minimum and maximum values to ensure that the input from this channel is negligible. The parameters listed below are derived from a sample 2 channel image.

```

inChNumberA = “00” #Red - YFP
inChNumberB = “03” #Green - GFP
inChNumberC = “00” #Blue – YFP (Placeholder channel)

```

```

sfA = 1.0
sfB = 1.0
sfC = 1.0

```

```

minA = 250
maxA = 800
minB = 200
maxB = 1500
minC = 10000
maxC = 10010

```

```

outChNumberH = “67”
outChNumberS = “66”

```

outChNumberB = “65”

2.4 Surface creation using the Brightness channel in Imaris and statistics extraction.

40. Start **Imaris**, and open the **original** TIFF files from
C:\Users\Lenovo\Desktop\Data\TIFF\P01_Renamed TIFFs.

► Troubleshooting

In later versions of **Imaris** with **Arena**, you may need to toggle to **Surfaces** mode first before you can access the **Open** button.

41. Press **Ctrl + I** to open the **Image Properties** panel of **Imaris**. Key in the following parameters for the voxel sizes, and press **OK** when done:

X: 0.155

Y: 0.155

Z: 1

▲ Critical step

It is important to key in the right values for the voxel size. If not, it will be necessary to regenerate the surface again in order to get the correct length scale.

42. If not already open, press **Ctrl + D** to toggle the **Display Adjustment** panel of **Imaris**. Set the following minimum and maximum values for the channels:

Channel 1: Min: 150 Max: 3500

Channel 2: Min: 220 Max: 3000

Channel 3: Min: 50 Max: 3500

□ Explanatory note

In the Base Color LUT, pixels with values smaller than or equal to the **Min** value is set as black (i.e. invisible), and pixels with values greater than or equal to the **Max** value is rendered in the maximum possible brightness (i.e. saturation – brighter pixels look no different from dimmer pixels beyond this point). By default, pixels with values in between are linearly interpolated and rendered with the corresponding hue intensity if **Gamma** is set as 1.0.

43. Press **Ctrl + I** again. Go to **Channels** panel, and in the **Base Color** tab, set the channel colors accordingly, and press **OK** when done:


Channel 1: Red: 0.000 Green: 0.000 Blue: 1.000 (i.e. pure blue)

Channel 2: Red: 0.000 Green: 1.000 Blue: 0.000 (i.e. pure green)

Channel 3: Red: 1.000 Green: 0.000 Blue: 0.000 (i.e. pure red)

□ Explanatory note

This sets the individual grayscale channels to utilize the base blue, green and red LUTs respectively. Thus, **Channel 1** appears blue, **Channel 2** green, and **Channel 3** red. You can also directly access the **Base Color** tab by clicking on the channel name in the **Display Adjustment** tab.



Channel color representation selection window of Imaris


□ Explanatory note

Take note of how **Imaris** represents colors using the **HSB model** in their color picker tool. If you key in a combination of values in the RGB model on the left, the specific color will be represented by the location of the small white square (in the hexagon) and rectangle (in the linear bar below). The six cardinal points of the hexagon represent the 'pure' hues (red, magenta, blue, cyan, green and yellow), and the angle of the white square relative to the origin indicates the corresponding **Hue** of the RGB combination. The relative distance from the origin (which is pure white) indicates its **Saturation** (i.e. how 'diluted' by white the RGB combination is). Finally, its position on the linear bar below indicates its **Brightness** (i.e. how 'dim' the color appears – note how the color bar goes to pure black when brightness is set to zero no matter what hue you set in the hexagon). If you are still confused by the HSB model, experiment with the color picker to get a better understanding of how it works.

44. Go to **Edit** > **Add** channels, and add the transformed **Hue**, **Saturation** and **Brightness channels** from *C:\Users\Lenovo\Desktop\Data\TIFF\P01_Renamed TIFs\P03_HSB*. The Brightness channel now appears as **Channel 4**, the Saturation channel is **Channel 5** and the Hue channel is **Channel 6**.
45. Set **Channel 4** as white using the similar procedure from **step 62** (**Red: 1.000** **Green: 1.000** **Blue: 1.000**). On the **Display Adjustment** tab, hide the rest of the


channels except for **Channel 4**, by unchecking the small check boxes beside the channel names.

46. In the **Surpass** tree, click **Add new surface** and use the following settings and go to the next step:




Surface creation step 1


47. Use the following settings and go to the next step:



48. Use the following settings and go to the next step:



49. Set the **Quality** value as **381**, and the following screen should appear:




Surface creation step 4

50. Go to the next step, and set the **Filter by voxel** value as **10**.
51. Click **Finish**.
52. To render the surfaces in their native hue, select **Surfaces 1** in the **Surpass** tree. Look for the **Color** tab, then go to **Color Type** box and select **Statistics Coded**. Under the **Statistics Type** drop-down option box, select **Intensity Median Ch=6**. Under **Colormap**, select the **HSB** look up table. Under **Colormap Range** set the minimum (**Min:**) as **0** and the maximum (**Max:**) as **40 000**.

□ Explanatory note

The surfaces are now rendered using the custom LUT that we earlier provided (**HSB.pal**). When fully stretched from **0** to **40,000**, the LUT maps the Hue channel (**Channel 6**) information to its expected color. The hue of the surfaces should reflect the actual colors of the cells when the Blue, Green and Red channels are simultaneously turned on (i.e. **Channels 1, 2 and 3**).



Surfaces rendered in their native hue


53. Save the **Imaris** file as **Testdata.ims** in **C:\Users\Lenovo\Desktop\Data** (click **File** > **Save As...**).

► Troubleshooting

In later versions of **Imaris** with **Arena**, you may need to click **File** > **Export** instead.

2.5 Saving out Surface Statistics

54. Go to **File** > **Preferences** > **Statistics**. Select the **Surfaces** tab and click on the checkbox until a tick is seen, as shown below.



Selection of types of statistics to include for Surfaces

□ Explanatory note

This step selects what kind of statistics is included to be calculated by Imaris. When the checkbox displays a tick, it means that all possible statistics will be included. If you wish, you can modify this by clicking on the small triangle on the left, and choosing only the relevant statistics you would like Imaris to report.

55. Select the surface (Surfaces 1), go to the Statistics tab, and select Detailed. Click on the button for **Export All Statistics to File** and save the statistics in a comma separated values (CSV) format into the following folder: *C:\Users\Lenovo\Desktop\Data* with the name “Testdata”. A subfolder *Testdata_Statistics* will now be generated in the selected folder containing multiple CSV files for each statistic. Close Imaris.

▲ Critical step

Avoid using any underscores in the name used for saving the statistics (excluding the ones automatically generated by Imaris itself) to prevent errors from occurring downstream when aggregating the statistics into a single CSV file.

2.6 Aggregate statistics into single CSV file and Gating in FlowJo

56. Run the script labelled *P04_Aggregate Imaris statistics output.py*. The following window should appear:



Graphical User Interface for P04

Set the following, save and close the file:

`prefix = "Testdata"`

□ Explanatory note

The variable `prefix` allows the user to customise the label for the file.

57. Click **Select directory**, and select the filepath `C:\Users\Lenovo\Desktop\Data\Testdata_statistics\` which contains the exported CSV files. Click **Done**. The script will consolidate the statistics from the individual files into a single CSV file.

```

C:\WINDOWS\system32\cmd.exe - python.bat "C:\Users\Lenovo\Desktop\Software\HSC\Scripts\P05_Aggregate Imaris statistics output_2017 12 06.py"
Reading Testdata_Intensity_Sum_Ch=4.csv...
S: Read in Testdata_Intensity_Sum_Ch=4.csv successfully... [Single]
Reading Testdata_Intensity_Sum_Ch=5.csv...
S: Read in Testdata_Intensity_Sum_Ch=5.csv successfully... [Single]
Reading Testdata_Intensity_Sum_Ch=6.csv...
S: Read in Testdata_Intensity_Sum_Ch=6.csv successfully... [Single]
Reading Testdata_Number_of_Disconnected_Components.csv...
S: Read in Testdata_Number_of_Disconnected_Components.csv successfully... [Single]
Reading Testdata_Number_of_Triangles.csv...
S: Read in Testdata_Number_of_Triangles.csv successfully... [Single]
Reading Testdata_Number_of_Vertices.csv...
S: Read in Testdata_Number_of_Vertices.csv successfully... [Single]
Reading Testdata_Number_of_Voxels.csv...
S: Read in Testdata_Number_of_Voxels.csv successfully... [Single]
Reading Testdata_Overall.csv...
Testdata_Overall.csv does not have Category, thus not included...
Reading Testdata_Position.csv...
Reading Testdata_Sphericity.csv...
S: Read in Testdata_Sphericity.csv successfully... [Single]
Reading Testdata_Time.csv...
S: Read in Testdata_Time.csv successfully... [Single]
Reading Testdata_Time_Index.csv...
S: Read in Testdata_Time_Index.csv successfully... [Single]
Reading Testdata_Volume.csv...
S: Read in Testdata_Volume.csv successfully... [Single]
C:\Users\Lenovo\Desktop\Data\Testdata_Statistics\P05_Aggregated statistics\Testdata_Objects.csv was created...
C:\Users\Lenovo\Desktop\Data\Testdata_Statistics\P05_Aggregated statistics\Testdata_Tracks.csv was created...
...Done. Press 'Enter' to exit.

```


Command window after P05 has completed running

58. A new subfolder, `P05_Aggregated_statistics`, containing 2 CSV files, `Testdata_Objects.csv` and `Testdata_Tracks.csv` will be created in the folder `C:\Users\Lenovo\Desktop\Data\Testdata_statistics\`.
59. Open **FlowJo**, and drag `Testdata_Objects.csv` into an empty workspace. It will automatically generate `Testdata_Objects.fcs` in the same folder as the CSV file, and load in the file.

▲ Critical step

FlowJo will automatically ignore a parameter where majority of the values are zero. If there are missing parameters in the drop-down menu, note down the missing parameter and check the values for the parameter in the original CSV file. Also note that for parameters whereby the values range from 0.0 to 1.0, such as Sphericity, the script P05 automatically multiplies the values by 10,000, such that the values now range from 0 to 10,000 instead. This is required for compatibility with FlowJo.

60. The **Testdata_Objects.fcs** sample should show up in the **FlowJo** workspace. Double-click on the sample filename, and the following window should show up.




FlowJo window of the **Testdata_Objects.fcs** sample

61. Set the X-axis as **Intensity Median Ch=5**, and the Y-axis as **Histogram**. Click the “**T**” button and select **Customise Axis** from the drop-down menu.


□ Explanatory note

Intensity Median Ch=5 is the intensity median of the saturation channel. In general, the median value is a more accurate reflection of the actual value as compared to the mean value as there may be some pixels with a high value, which may skew the actual mean value.




Setting display type and range of the axes scales in FlowJo

62. The **Scale** should be set as **Linear**, the **Min** value as **0**, and the **Max** value as **6500**. Click **Apply**.
63. Use the **Range** gating button (beside the arrow icon on the top left corner) to set the minimum and maximum limits of the gate as shown. You should now obtain a gate for the triple positive cells. Label the population as **“White”**.




Histogram gating of the triple positive ('White') population

64. In order to identify the spatial localization of the gated cells, go to the workspace and double click on the **White** population.




Selecting the 'White' population in the FlowJo workspace

The FlowJo window should open, set the X-axis as Pos X and the Y-axis as Pos Y. As there are very few surfaces, it is difficult to see the individual dots. To overcome this, click Display and select the option for Draw Large Dots.




Selection of 'Draw Large Dots' option in FlowJo

65. To export the graph, right click the window and select **Copy to Layout Editor**. The **Layout Editor** should pop up. The following image should be obtained:



Exported spatial localization from Layout Editor in FlowJo

66. The population can now be exported for backgating. In the **FlowJo** workspace, right-click the **White** population and select **Export/Concatenate Populations**. The following window should pop up.



FlowJo export of populations

Use the settings as shown above. Take note that the **Format** used should be **CSV-Scale values** so that the ID of the surface objects are not modified by **FlowJo**. Set the **Destination** as **C:\Users\Lenovo\Desktop\Data\Population**. Click **Export**.

▲ Critical step


Do not choose **CSV-Channel values**, as the numbers within each channel will be scaled by **FlowJo** and thus become changed. Theoretically, this scaling function should only happen when the export format of **CSV-Scale values** is chosen, but the scaling function has been swapped in this version of **FlowJo**. This is a bug that has yet to be rectified by the software manufacturer, but may be corrected in future versions.

2.7 Backgating in Imaris (Version 8.4.2)

67. Open the **Imaris** file **Testdata.ims** in **C:\Users\Lenovo\Desktop\Data**. Ensure that the mouse cursor has selected the surface (**Surfaces 1**) as shown below.


▲ Critical step

If the correct surface is not selected, there will be an error message during the running of the script since the script is acting on the wrong object.



Backgating in Imaris Step 1

68. Go to **Image Processing** > **XTIsolateSurfaceswithCsvList** and start the script. The following dialog box will pop up.




Backgating in Imaris Step 2

69. Click **OK**. Go to the filepath *C:\Users\Lenovo\Desktop\Data\Population* and select **export_Testdata_Objects_White.csv**, then click **Open**.
70. A new folder labelled **Surfaces 1_subset** should appear in the **Surpass tree**. Select the surface called **Surfaces 1_Select** which was created in this new folder. The color of the surfaces can then be changed by going to the **Color** tab, looking for the **Color Type** box, and selecting **Base**.
71. To check if triple positive cells have been correctly gated, set **Channel 5** (saturation channel) as white using the similar procedure from **step 62** (**Red: 1.000** **Green: 1.000** **Blue: 1.000**). On the **Display Adjustment** tab, overlay **Channel 5** together with

the original three channels by checking **Channels 1, 2, 3** and **Channel 5** using the small check boxes beside the channel names. Set the following minimum and maximum values for **Channel 5**: **Min: 0, Max: 4250**.

72. To ensure that the signal from the original 3 channels can be viewed, uncheck **Surfaces 1**. The screen should show something similar to what is shown below. To turn the backgated surfaces on and off, check and uncheck **Surfaces 1_Select**.



Backgated surfaces (right) compared to Saturation channel overlaid on the original 3 channels (left). These cells were originally very difficult to identify by eye alone.

▲ Critical step

XTIsolateSurfaceswithCsvList works with **Imaris 8.4**, and as of the date of this tutorial creation, is not yet able to work with **Imaris 9.0** or above. We were informed by the company representative that the developers are still working on getting Python scripts to work with surfaces in **Imaris 9.0**.

► Troubleshooting

If the population contains only 1 surface, the script will not be able to run. However, in this case, there should not be any point in running the script at all, but if absolutely required (e.g. for consistency), the user can perform a workaround by manually duplicating the surface of interest first.

► Troubleshooting


If an error occurs during the backgating, the user should double-check that the exported surface IDs from **FlowJo** match the original IDs of the surface. If the exported surface IDs from **FlowJo** are different, it could be due to a bug in **FlowJo**. You can try closing all the programs completely and repeating **Section 2.7** again.

3 Spectral Compensation (optional)

3.1 Compensation using FlowJo

In the case where there is spillover, spectral compensation is required. For the purposes of this tutorial, we will be using a test data set with a mixture of YFP- and GFP-containing mostly single-positive cells. Since we can unambiguously identify our single-positive cells using the data set itself, we do not need any ‘single stain’ controls for this example, but you may need to prepare them for your own dataset. You should acquire the ‘single stain’ controls on the same day with the exact same settings as your actual data.

73. Copy the folder called *C:\Users\Lenovo\Desktop\Software\HSC\Data\Compensation* to the filepath *C:\Users\Lenovo\Desktop\Data\Compensation*.
74. Drag the file *Compensation_Objects.csv* from the *C:\Users\Lenovo\Desktop\Data\Compensation* folder into *FlowJo*. The program should automatically convert the CSV file into a FCS format.
75. Plot *Parameter_6: Intensity_Center_Ch_01* (YFP) for the X-axis, and *Parameter_7: Intensity_Center_Ch_02* (GFP) as the Y-axis. Click the **T** button for the X-axis. The scale should be set to **Linear**. Select **Customise Axis...** and set the minimum value as **50**, maximum value as **2000**. Do the same for the Y axis, but set the minimum value as **0**, maximum value as **2500**. See **Step 81-82** to recap if you have forgotten how to do this.
76. Select **Display** and click **Draw Large Dots**. Click on the **Polygon** gate button (four buttons away from the arrow icon on the top left corner) and gate the YFP⁺ and GFP⁺ cells as shown below:



Individually gated populations

77. Create a new subfolder *Compensation_Statistics* under *C:\Users\Lenovo\Desktop\Data\Compensation*. Go to the **FlowJo** workspace, select both the YFP and GFP populations and right click. Select **Export/Concatenate Populations**. Export both the populations as **CSV-scale values** in this folder. Click **Export**. The filenames that should appear are *export_Compensation_Objects_YFP.csv* and *export_Compensation_Objects_GFP.csv*.

► Troubleshooting

In older versions of **FlowJo**, you may need to export the two populations individually.

78. Run **Microsoft Excel**. Go to **File > Open** and point to the file *export_Compensation_Objects_YFP.csv*. Keep the **FlowJo** workspace open.

► Troubleshooting

Depending on your version of **Microsoft Excel**, you may need to open the file with a different procedure. Also, you may need to change the file type to the **All files (*.*)** option in order to see the CSV files in the folder.


79. Select **Column R** (*Intensity_Center_Ch_01*) (YFP) and **Column S** (*Intensity_Center_Ch_02*) (GFP), and press **Ctrl + C** to copy the columns. Create a

new worksheet, select **Cell A1** (uppermost left corner), then press **Ctrl + V** to paste the data onto the new worksheet.

► Troubleshooting


Depending on your version of **FlowJo**, the exact columns for **Intensity_Center_Ch_01** and **_02** may not be the same. Adjust accordingly if this is the case.

80. Select **Columns A and B**. Click **Insert > Scatter** to plot **Intensity_Center_Ch_02** (Y-axis) against **Intensity_Center_Ch_01** (X-axis). Select the **Scatter** (without lines) type of graph as shown in the screenshot below.



Insert scatter plot without lines

81. Left-click the graph to select the datapoints, right-click and select **Add Trendline...**. The **Trendline Options** window should pop up. Check the option box for **Display Equation on chart**.



Calculating the coefficient for YFP spillover into GFP channel

Based on the gradient of this trendline, we can see that the percentage spillover from the YFP cells into the GFP channel is approximately **35.19%**.

▲ Critical step

Your aim is to calculate the slope as accurately as possible, and you should remove any outliers that you can spot. In our example above, the data point at (1684.39, 552.93) is clearly out of place, and when removed would modify the slope to 38.52%. For the sake of simplicity, we did not include this step in the tutorial, but you should be aware of this.

▲ Critical step

Depending on the acquisition parameters and the excitation/emission dye spectra, note that there may not necessarily be any spillover. If the graph does not have an obvious correlation, or if you get a negative gradient, the true gradient is likely to be zero. You can turn on the **Display R-squared value on chart** option in the **Trendline Options** window to show the R^2 statistic, which ranges from 0.0 to 1.0. A very low value indicates weak correlation, and may indicate no spillover in reality.

▲ Critical step

To check for spillover from YFP cells into the GFP channel, it is critical that **Intensity_Center_Ch_01 (YFP cells)** is plotted as the X-axis, and that **Intensity_Center_Ch_02 (detected spillover signal in GFP channel)** is plotted as the Y-axis. Plotting the values in the wrong manner will lead to the wrong coefficients.

□ Explanatory note

The calculation for the spillover coefficient is defined by the equation $y = mx + c$, where y represents the signal detected in the channel of interest and x is the signal from the original spilling channel. The spillover coefficient is the gradient m , which is the percentage of the signal from the original channel that has spilt into the detector of the other channel of interest. In this case, since the source cell population is YFP+, **Channel 01** is the spilling channel, and thus always represents x . If you have more than 2 channels, you would modify this step to plot various channels as y , but always keep **Channel 01** as x (e.g. **Ch_02**, **Ch_03**, **Ch_04**... vs **Ch_01**). Take note that the compound offset value c is inconsequential for subsequent linear algebra calculations, and only the slope m is critical.

□ Explanatory note


We strongly recommend using the **Intensity_Center_Ch_XX** parameters for deriving the spillover coefficients, instead of other aggregate parameters such as intensity mean or median values. This parameter type refers to the channel intensities of the same pixel at the center of the surface, and thus most accurately reflect the spillover profile of the images on a per pixel basis.

82. Go to **File > Open** and point to the file **export_Compensation_Objects_GFP.csv**.
83. Select **Column R (Intensity_Center_Ch_01)** (YFP) and **Column S (Intensity_Center_Ch_02)** (GFP), and press **Ctrl + C** to copy the columns. Create a new worksheet, select **Cell A1** (uppermost left corner), then press **Ctrl + V** to paste the data onto the new worksheet.
84. Now, however, we need to ensure that **Intensity_Center_Ch_02 (GFP cells)** is recognized as the X-axis in this case, and **Intensity_Center_Ch_01 (detected spillover signal in YFP channel)** is recognized as the Y-axis. To do this, select **Column A** (YFP signal), press **Ctrl + X** to cut the column, then select **Column C**, and press **Ctrl + V** to paste.

□ Explanatory note

By default, Microsoft Excel always recognizes the left-most column as the X-axis.


85. Plot **Intensity_Center_Ch_01** (Y-axis) against **Intensity_Center_Ch_02** (X-axis) in the manner described in **Steps 100 & 101**. This time, you should select **Columns B and C** in this case. A graph similar to what is shown below should be seen:



Calculating the coefficient for GFP spillover into YFP channel

The amount of GFP signal spillover to YFP is given by the gradient, which is the value **0.2276**. In percentage terms, it is **22.76%**.

86. In the **FlowJo** workspace, select **Tools** and **Edit Compensation Matrix**. The following window should pop up. Select the **+** button at the bottom left hand corner of the window and **New Identity Matrix** to add a new matrix. A window **Choose Selected Parameters** will pop up, select **Parameter_6: Intensity_Center_Ch_01** and **Parameter_7: Intensity_Center_Ch_02**. Hold the **Shift** key to select multiple parameters. Click **Choose Selected Parameters**.



Add new Compensation Matrix

87. Key in the spillover coefficients into the compensation matrix in the format as shown below:


	Parameter_6 :: Intensity_Center_C...	Parameter_7 :: Intensity_Center_C...
Parameter_6 :: Intensity_Center_C...	100	35.19
Parameter_7 :: Intensity_Center_C...	22.78	100

Compensation Matrix in FlowJo

□Explanatory note


The **FlowJo Compensation Matrix** is formatted in a manner where the **dyes** are represented by the **rows** and the **detectors** are represented by the **columns**. Hence, the spillover coefficient (35.19%) for YFP cells (**Parameter_6**, first row) into the GFP detector (**Parameter_7**, second column) is keyed in the first row, and vice versa for GFP cells (**Parameter_7**, second row) into the YFP detector (**Parameter_6**, first column). By definition, the amount of signal generated by YFP cells (first row) detected by the YFP detector (first column) is 100%, and the same holds true for GFP cells (second row) into the GFP detector (second column). The values in the diagonal are thus always fixed at 100.

88. Click the **[M]** button, hold, drag and drop it over the **Compensation_Objects.fcs** in the **FlowJo** workspace. A colored grid should now appear beside the **Compensation_Objects.fcs** file (indicated by the red arrow in the figure below). This indicates that the compensation matrix has been applied to the sample.



Applying the compensation matrix to *Compensation_Objects.fcs*. When a compensation matrix has been applied successfully to a sample, a colored grid should show up beside it, as indicated by the red arrow here.

89. In order to check that the GFP and YFP cells have been accurately corrected, go to the *Compensation_Objects.fcs* plot window, select **Display > Show Uncompensated Parameters**. Select the X-axis as **Parameter_6: Intensity_Center_Ch_01** and the Y axis as **Parameter_7: Intensity_Center_Ch_02**. Right click and select **Copy to Layout Editor**. Next, select the X-axis as **Comp-Parameter_6: Intensity_Center_Ch_01** and the Y axis as **Comp-Parameter_7: Intensity_Center_Ch_02**. Right click and select **Copy to Layout Editor**. An image like the one below should appear:



Comparison of dotplots of uncompensated and compensated channels

Note the perpendicular angle between the YFP and GFP dots, indicating a successful compensation.

If you do not need to correct the image, but just the histocytometry plots, then you can stop here. If the images themselves need to be corrected, then you should proceed with the steps in the next subsection.

3.2 Correction of original images

90. Run the script **S01_Channel unmixing.py**. The following window should pop up:



Graphical user interface of script S01

91. Key in the following values for the input channel numbers:

inChNumbers = 00 01

□ Explanatory note

Any number of channels can be corrected by this script. In this case, channel “00” is YFP, channel “01” is GFP.

92. Key in the following values for the compensation matrix, then save and close **Notepad**:

cMatrix = [[1.0, 0.3519], [0.2276, 1.0]]

□ Explanatory note

The matrix should be typed in the **FlowJo** format where the **row** headers are **Dyes**, **column** headers are **Detectors**. For this script, use ratios instead of percentages (e.g. use 1.0 instead of 100%).

93. Run the script, and point to the original uncorrected channels in **C:\Users\Lenovo\Desktop\Data\Compensation\CompensationTIFF** to generate the corrected channels in the subfolder **S01_Unmixed**.
94. Open **C:\Users\Lenovo\Desktop\Data\Compensation\Compensation.ims** with **Imaris**, which already has the surfaces of the cells. Ensure that **Surfaces 1** is unchecked. To load in the corrected channels into the file, go to **Edit**, and select **Add Channels....** Point to **C:\Users\Lenovo\Desktop\Data\Compensation\CompensationTIFF\S01_Unmixed** and click **Open**. Both the uncompensated and compensated channels will now be present. If the **Display Adjustment** tab is not open, click **Ctrl + D**. If the **Display Adjustment** tab is open, use the following settings to visualise the uncompensated and compensated channels using the following colours:


Channel 1 (uncompensated YFP)	Min: 250	Max: 950 - red
Channel 2 (uncompensated GFP)	Min: 250	Max: 1250 - green
Channel 3 (compensated YFP)	Min: 210	Max: 910 -yellow
Channel 4 (compensated GFP)	Min: 610	Max: 1610 -blue

Click **Channel 3** in the **Display Adjustment** tab and under the **Base Color** tab, change the settings for **Channel 3** and **Channel 4** to:

Channel 3	Red: 1.000	Green: 1.000	Blue: 0.000 (yellow)
Channel 4	Red: 0.000	Green: 0.000	Blue: 1.000 (blue)

Toggle the various combinations of channels to see the effects of the channel unmixing on the original data.

You should be able to get the images illustrated below:



Comparison of compensated and uncompensated channels.

The GFP signal from the neutrophils that spilt into the YFP channel have been removed and the YFP signal from the dendritic cells that spilt into the GFP channel have been removed in the compensated channels. Now that the channels have been appropriately compensated, the user can then proceed with any downstream steps that they need to do for their own purposes.